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1. DESCRIPTION OF THE IMPLANT
Textile structure (knitted) in polypropylene (PP)

2. PERFORMANCES
Uro-gynaecological support.

3. OPERATING TECHNIQUE
The surgeon decides what operating technique to use according to his or her experience in this field.

4. INDICATIONS
· Stress urinary incontinence

5. CONTRA-INDICATIONS
· Acute or chronic, local or systemic infections.
· Severe muscular, neurological, metabolic or vascular deficiency at the affected site .
· Any accompanying disorder liable to affect the function of the implant
· Pregnancy
· Growing patient
· Allergy or reaction to PP

EN

6.PRECAUTIONS FOR USE
DO NOT REUSE-DO NOT RESTERILISE
·The device must be stored in the protective case provided for each model. The case must never be opened before the operation starts.
·The device must not be used if the primary packaging is damaged. Sterility cannot be ensured if the case is damaged. Check it carefully before use.
·Examine each implant before use to detect any flaws (cracks, distortion, etc.) and reject it if necessary.
·The device must be implanted using the right ancillary provided with the implant.
- Remove the ancillary respecting its path of insertion and its curvature.

7. POSSIBLE SIDE-EFFECTS
· Intraoperative complications: Wounds or damage to surrounding organs during surgery (bladder, urethra, vagina, rectum, blood vessels), particularly during retropubic surgery.
· Specific early complications (< 6 weeks): Perineal or pelvic haemorrhage/hematoma; Acute urinary retention (particularly in the case of transobturator surgery); Acute postoperative pain; Postoperative infections.
· Specific late complications (> 6 weeks): Vaginal, bladder or urethral prosthetic exposure to be looked for in case of: Chronic pain, Bladder hyperactivity, Recurrent urinary tract infections, Haematuria, Dyspareunia (pain during intercourse),

Abnormal vaginal bleeding or discharge; Bladder hyperactivity in case of obstruction or urinary infection; Late dysuria, particularly after retropubic surgery; Chronic pain at the root of the thigh, particularly after transobturator surgery; or at the level of the passage of the sling.

STORAGE/EXPIRY DATE
Do not use the device after the last day of the month of expiry.

TRACEABILITY
A traceability label identifying the batch and serial numbers of the device is available inside the box. This label is intended to be attached to the patient's medical records.

IMPORTANT: THE DEVICE MUST ONLY BE USED BY AN ADEQUATELY-TRAINED AND EXPERIENCED UROGYNÆCOLOGICAL SURGEON.
THE SURGEON MUST INFORM THE PATIENT OF THE CONTRA-INDICATIONS AND POSSIBLE SIDE EFFECTS.
THE SURGEON MUST INFORM THE PATIENT OF THE POST-OPERATIVE PROCEDURE TO BE FOLLOWED AND OF THE RISK OF FAILURE IF IT IS NOT FOLLOWED

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