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## ENGLISH VERSION

### 4. OBJECT

This technical note explains how to reprocess and keep the leak test socket in the bowl clean. It also describes the daily maintenance actions to be carried out. This will prevent any source of contamination and ensure optimal hygiene.

**Reminder:** Failure to comply with the maintenance instructions may invalidate the manufacturer's warranty and can adversely affect the quality of the device reprocessing.

### 5. DAILY MAINTENANCE

Daily maintenance is required to ensure optimal performance of Soluscope washer disinfectors. These actions must also be performed after each maintenance intervention on the machine.

- Wipe the outside of the washer disinfectant with a pH-neutral cleaning and disinfecting wipe.
- Clean the washer disinfectant plate, the interior edge of the bowl (first 5 centimeters from the edge) and the lid, using a surface disinfectant spray (e.g. Aniospray Surf 29 by Laboratoires Anios).
- Reprocess the leak test socket (*figure 1*) in the bowl (see reprocessing of the bowl leak test procedure paragraph 3).
- Clean the lid seal using a surface disinfectant spray.
- Check that there is no water in the leak test connector. Use a clean cloth soaked in 70% ethyl alcohol to disinfect the SL-OTE connector tip (*figure 2*). Make sure the connector is completely dry before plugging it in.
- Launch self-disinfection of the washer-disinfectant (cycle 6).
- Dry and check the SL-OTE connectors (*figure 2*), then disinfect them using a surface disinfectant spray. Make sure there are no cracks or kinking marks. Store the connectors in a clean and dry location. Lubricate the seal of leak test connectors with silicone oil (Lubrispray P recommended).

### 6. METHOD

#### 6.1. PREPARATION

Wear gloves to protect your hands.

Make sure the 70% ethyl alcohol is in a ready-to-use spray bottle (not diluted in-house).

Disconnect all leak test connectors.

**Important reminder:** Using a water-based disinfectant as an alternative is strictly prohibited. It may damage the endoscope.

## 6.2. INITIAL CLEANING

If the leak test socket is visibly soiled, be sure to clean the inside of the socket with a swab or brush to remove any dirt or debris.

*Figure 1: Leak test socket in the bowl*

## 6.3. SPRAYING

Spray (2 sprays) 70% ethyl alcohol directly onto the leak test socket from approximately 10cm. Be sure to cover all internal surfaces of the leak test socket.

## 6.4. CONTACT TIME

Leave for 5 min to ensure effective action. Depending on the disinfectant used, please follow the contact time recommended by the manufacturer.

## 6.5. BLOWING

### a) For software versions up to v1.8.4

Do not connect anything on leak test socket in the bowl. Close the lid and launch cycle 11. After 2 minutes of blowing, an alarm 1 should appear. Return to main view.


### b) For software versions v1.9.0 and above

Use cycle 10 to blow off any excess alcohol and wait for the end of cycle. At this step, wipe the lid seal.

## 6.6. WIPING

Use a clean cloth soaked in 70% ethyl alcohol to reprocess the SL-OTE plug tip.

Make sure the plug is completely dry before plugging it in.

*Figure 2: SL-OTE area to be reprocessed*

### **6.7. PRECAUTIONS**

Avoid using too much alcohol to avoid flooding the leak test socket (Do not spray on contact with leak test socket).

### **6.8. WDE SELF-DISINFECTION**

Connect SL-OTE plug to leak test socket on bowl and launch cycle 6.